



OPERATION THEATER TABLES (O.T. TABLES)

A). C –Arm Provision Hydraulic O.T. Table



Technical Specification

Minimum Height	:	75 cm
Maximum Height	:	105 cm
Length of Table Top	:	195 cm to 215 cm
Width of Table Top	:	50 cm
Lateral Tilt	:	$20^{\circ} \pm 2^{\circ}$
Trendelenburg	:	$30^{\circ} \pm 2^{\circ}$
Chair Position	:	$-22^{\circ} / +90^{\circ}$
Head Section	:	$+90^{\circ} / -90^{\circ}$
Leg Position	:	$+90^{\circ} / -45^{\circ}$
Kidney Elevation	:	15 cm

09898622773

www.sunmedicares.com

contact@sunmedicares.com

B - 5 "Kailas", B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



Salient Feature :

Compatible with C arm image intensifier used for orthopedics surgery.

Can also be used for general surgeries like Cardiac, Urology, Endoscopy. The table top is made of radio-translucent material.

Head and Leg sections are Interchangeable.

Kidney Bridge detachable or inbuilt as per customer requirement.

Base and Center column are made of Stainless Steel or Fiber Glass for easy cleaning and hygienic.

Table frame, side rails and accessories are made up of stainless steel.

Accessories of O.T. Table:

Anesthetic Frame (1 Pcs)

Padded Side Support (1 Pair)

Intravenous Armboard (1 Pair)

Padded Shoulder Support (1 Pair)

Knee Crutches (1 Pair)

One Mattress Set - Thickness 50 mm

Note:

Option Provided - Up-Down Movement Electrically, Rest position through manual hand operated crank system. Also Option for Longitudinal Top Side : 20 Cm



09898622773



www.sunmedicares.com



contact@sunmedicares.com



B - 5 "Kailas" , B /h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



VARIOUS POSTION OF O.T.TABLE



Trendelenburg



Reverse Trendelenburg



Chair Position



Flex Position



Reflex Position



Lithotomy Position



Kidney Position



Lateral Tilt(Both)



Back Position



B). ELECTRIC O.T. TABLE WITH MANUAL OVERRIDE **(Electric operated, Battery Operated, Manual Operated)**



Technical Specification

Minimum Height	:	75 cm
Maximum Height	:	105 cm
Length Of table Top	:	215 cm
Width Of Table Top	:	50 cm
Lateral Tilt	:	+20° / -20°
Trendelenburg	:	+25° / -25°
Back Section	:	+90° / -40°
Head & Leg Section	:	+90° / -90°
Kidney Elevator	:	15 cm manual operated

09898622773

www.sunmedicares.com

contact@sunmedicares.com

B - 5 "Kailas" , B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



Salient Feature :

100% Life Saving due to voltage Electricity Supply.
Battery Back-up.
800 lb. weight capacity.
Fully Table Function, Fully Articulation Capability.
Removable Leg Section.
Full Body Imaging without special accessories.
Durable one-touch pendant control.
Casters lockable with single action brake.
All stainless Steel Accessories.
Electrically Controlled Height, Backrest, Trendelenberg / Reverse, Lateral Tilt.
Radiolucent table top.
Manual over ride by rotating crank on base.
Head, Backrest, Trunk and Leg 4 – Section table top.
5 – Section Seamless molded mattress.
Support on a heavy base column covered with stainless steel sheet.

(Option available - Top Side Both Head & Leg Side: 30 cm manual operated.)

Accessories of O.T. Table:

Anesthetic frame (1 Pcs)
Padded Side Support (1 Pair)
Intravenous Arm board (1 Pair)
Padded Shoulder Support (1 pair)
Knee Crutches (1 Pair)
One Matters Set Thickness – 50 mm

09898622773

www.sunmedicares.com

contact@sunmedicares.com

B - 5 "Kailas" , B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



VARIOUS POSTION OF O.T.TABLE



Lateral Tilt (Both)



Trendelenburg



Reverse Trendelenburg



Chair Position



Flex Position



Reflex Position



Kidney Position



C). OPHTHALMIC O.T. TABLE (Electrically , Battery & Manually Operated)



Technical Specification

Length of Table Top	: 2100 mm
Width of Table Top	: 500 mm
Minimum Height	: 650 mm (10%±) Motorized
Maximum Height	: 950 mm (10%±) Motorized
Back Section Up	: 90° Motorized
Back Section Down	: 30° Motorized
Trendelenburg	: 25° Motorized
Reverse Trendelenburg	: 25°

Salient Features:

100 % Life Saving Due to low Voltage Electricity Supply.
Offline Battery backup.
800 Lb. weight capacity.
Full Table Functions, Full Articulation Capability.
Durable one-touch pendant control.
Casters lockable with single action brake.
All stainless steel accessories.
Electrically Controlled Height, Trendelenburg / Reverse,
Lateral tilt, Head, back & leg section.
Radiolucent Table Top.
Support on a heavy base column covered with stainless steel sheet.



09898622773



www.sunmedicares.com



contact@sunmedicares.com



B - 5 "Kailas", B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



D). FULLY AUTOMATIC O.T TABLE (Electrically, Battery , Manual Operated)



Trend Rev



Reverse Trend Rev



Flex



Reflex



Tilt



Chair



E). Gynaec O.T. Table – Hydraulic Operated.



Technical Specification

Minimum Height	:	700 mm
Maximum Height	:	950 mm
Length of Table Top	:	1870 mm
Width of Table Top	:	610 mm
Lateral Tilt	:	20° ±2°
Trendelenburg	:	30° ±2°
Hydraulic Lift	:	10"
Back position	:	+ 90 ° / - 40°
Kidney Position	:	120 mm
Weighting Capacity	:	200 Kg

09898622773

www.sunmedicares.com

contact@sunmedicares.com

B - 5 "Kailas" , B /h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



Salient Feature

- : Single Action Hydraulic Lift.
- : Heavy Mild Steel Base with fiber Wheels.
- : Auto Floor Locking Devices.
- : The Detachable Leg Section Provide Easy Lithotomy Positioning.
- : The improved hydraulic System Provides Smooth Movment And the Speed of Lowering Permanently Set.
- Head and Leg Section interchangeable.
- Weight Bearing capacity 200 Kgs.

Accessories

Anesthetic Screen Frame with Clamp	: 1 pcs.
Padded Side Supports	: 1 pair
Intravenous Arm board	: 1 pair
Padded Shoulder Support	: 1 Pair
Knee Crutches	: 1 pair
Mattress Set Thickness 50 mm	: 1 set
Handle	: 1 Pair



09898622773



www.sunmedicares.com




contact@sunmedicares.com



B - 5 "Kailas" , B /h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



F). Gynaec ELECTRIC O.T. Table (Battery Back Up)



Technical Specification

Minimum Height	: 635 mm
Maximum Height	: 965 mm
Backboard fold upward	: 80°
Length of the table	: 1750 mm
Width of the table	: 600 mm

Accessories Of O.T. Table

Anesthetic Frame (1 Pcs)
Padded Side Support (1 Pair)
Intravenous Arm board (1 Pair)
Padded Shoulder Support (1 Pair)

Salient Feature

100% Life Saving Due to low voltage electricity supply
Battery Backup.
800 lb. weight capacity
Full Table functions, full Articulation capability.
Removable leg section.
Durable One touch pendant control.
Caster Lockable with single action brake.
All stainless steel accessories.
Electrically Controlled Height
Radiolucent Table Top.

09898622773

www.sunmedicares.com

contact@sunmedicares.com

B - 5 "Kailas", B /h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.