



A). VENTILATOR - Portable / Transportable / Ambulance Use / Battery Backup (Stand Model also Available) TURBINE BASE VENTILATOR.



(A TRUE LIFE SAVER IN THE BATTLE OF INTESIVE CARE)

TECHNICAL SPECIFICATION

AVAILABLE I.V. MODES	: VCV, PCV, VACV, VSIMV, PSIMV, PACV, CPAP
AVAILABLE N.I.V MODES	: CPAP, BiPAP, HFT
PEEP	: 0 TO 30 CM H ₂ O
RASPIRATORY RATE	: 1 TO 80 BPM
INSPIRATORY PRESSURE	: 5 TO 100 CM H ₂ O (HIGH PRESSURE LIMIT – LOW PRESSURE LIMIT)
TIDAL VOLUME	: 40 TO 2000 ml
INSPIRATORY TIME	: 0.1 TO 9.9 Sec.
I : E RATIO	: 1 : 9.9 - 2. 5:1
TRIGGER	: PRESSURE (0 TO -20 CM H ₂ O)
FIO ₂	: 21% TO 100%
DIMENSIONS (WX D X H)	: 300 mm x 260 mm x 130 mm
EX SENCE	: 50% TO 100%
APNEA BACK UP	: YES
DETACHABLE AND AUTO CLAVABLE	: YES
EXPIRATORY VALVE	
RECHARGABLE Li-ion BATTERY	: YES WITH UPTO 12 HOURS BACK UP



09898622773



www.sunmedicares.com



contact@sunmedicares.com



B - 5 "Kailas", B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



B). VENTILATOR - Screen 7" TFT (Stand Model also Available) - TURBINE BASED ICU VENTILATOR

**(A TRUE LIFE SAVER IN
THE BATTLE OF ICU CARE)**



TECHNICAL SPECIFICATION

NEBULIZER	: YES
DISPLAY	: 7" TFT CAPACITIVE
AVAILABLE I.V. MODES	: VCV, PCV, VACV, VSIMV, PSIMV, PACV, CPAP
AVAILABLE N.I.V MODES	: CPAP, BiPAP, HFT
WAVEFORMS	: VOLUME, PRESSURE, FLOW
PEEP	: 0 TO 30 CM H ₂ O
RASPIRATORY RATE	: 1 TO 80 BPM
INSPIRATORY PRESSURE	: 5 TO 100 CM H ₂ O (HIGH PRESSURE LIMIT – LOW PRESSURE LIMIT)
ENCODER	: YES
TIDAL VOLUME	: 40 TO 2000 ml
INSPIRATORY TIME	: 0.1 TO 9.9 Sec.
I : E RATIO	: 1 : 9.9 - 2.5 : 1
TRIGGER	: PRESSURE (0 TO -20 CM H ₂ O)
FIO₂	: 21% TO 100%
DIMENSIONS (WX D X H)	: 310 mm x 350 mm x 280 mm
EX SENCE	: 50% TO 100%
APNEA BACK UP	: YES
DETACHABLE AND AUTO CLAVABLE	: YES
EXPIRATORY VALVE	
RECHARGABLE Li-ion BATTERY	: YES WITH UPTO 12 HOURS BACK UP

09898622773

www.sunmedicares.com

contact@sunmedicares.com

B - 5 "Kailas", B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



C). VENTILATOR - Screen 10" TFT - TURBINE BASED ICU VENTILATOR

**(A TRUE LIFE SAVER IN
THE BATTLE OF ICU CARE)**

TECHNICAL SPECIFICATION

NEBULIZER	: YES
DISPLAY	: 10" TFT CAPACITIVE
AVAILABLE I.V. MODES	: VCV, PCV, VACV, VSIMV, PSIMV, PACV, CPAP
AVAILABLE N.I.V MODES	: CPAP ,BiPAP , HFT
WAVEFORMS	: VOLUME, PRESSURE, FLOW
PEEP	: 0 TO 30 CM H ₂ O
RASPIRATORY RATE	: 1 TO 80 BPM
INSPIRATORY PRESSURE	: 5 TO 100 CM H ₂ O (HIGH PRESSURE LIMIT – LOW PRESSURE LIMIT)
ENCODER	: YES
TIDAL VOLUME	: 40 TO 2000 ml
INSPIRATORY TIME	: 0.1 TO 9.9 Sec.
I : E RATIO	: 1 : 9.9 - 2. 5:1
TRIGGER	: PRESSURE (0 TO -20 CM H ₂ O)
FIO₂	: 21% TO 100%
DIMENSIONS (WX D X H)	: 310 mm x 350 mm x 280 mm
EX SENCE	: 50% TO 100%
APNEA BACK UP	: YES
DETACHABLE AND AUTO CLAVABLE	: YES
EXPIRATORY VALVE	
RECHARGABLE Li-ion BATTERY	: YES WITH UPTO 12 HOURS BACK UP

