



FOR TROLLEY WASH AREA



1 S.S. Unsterile Closed Transportation Collection Trolley

Size :- 950 x 750 x 1150 mm LxWxH



2 Trolley Wash System

Portable or Fix Wash System provided as per the criteria. In General Cart / Trolley Wash in CSSD can be done.



Manual:

- a). Water Pressure driven between 2 to 6 Bar up to 50° C.
- b). Dosing of Detergent should be Automatic.
- c). Wall mounting Kit.
- d). Cleaning Gun 15 mm
- e). Length hose with Stainless Steel ½" to ¾" connection.

Automatic:

- a). An Automatic mechanical cleaning equipment that help reduce microbial contamination and render items safe for staff to handle.

Using Both Techniques:

In the case the manual carts wash area is serving as a back-up to the automated cart washer.

3 Air Gun with one Single Nozzle and with Pipe assembly



4 Water spray Jet Gun with one Single Nozzle and with Pipe assembly






FOR WASH AND DECONTAMINATION AREA



1 S.S. Soiled Receiving Table
Size :- 1200 x 600 x 900 mm LxWxH



**2 S.S. Double/Triple Sink Wash Station with
2/3 Deep Sinks and with other Accessories**
Unit Size:- 2000 x 650 x 900 + 150 mm LxWxH + Flashback
Sink Size:- 600 x 450 x 250 mm LxWxDeep



**3 Air Gun with one Single
Nozzle and with Pipe assembly**



5 Imported Air / Water Spray Gun with 8 different Nozzle attachments and with Wall mounted Rack

The Medical Spray Guns are designed for use within soiled pre-cleaning areas in reprocessing units. Especially suitable for fine-mesh instruments. The powerful, multi-purpose spray-gun rinse uses water or compressed air to rinse your items. The spray-gun is easy to operate with one hand, enabling the operator to find the best angle for rinsing. The heat-insulated gun grip provides extra safety. Air Consumption 8-9 CFM, Air Pressure : 70-100 PSI. Forward Speed :Adjustable – Constant



6 Eye Wash Station :



7 S.S. Working Table

Size :- 1200 x 600 x 900 mm LxWxH



8 S.S. Free standing Open Storage Rack with 5 Heavy duty wire.
Shelves Size :- 1200 x 450 x 1800 mm LxWxH



9 S.S. Pass Box - Double Door with Manual Mechanical Interlocking and with Transparent Vision Panel .
External Size :- 600 x 650 x 600 mm LxWxH




FOR PACKING AND PREPARATION AREA:



1 S.S. Instrument Receiving Working Table

Size:- 1200 x 600 x 900 mm LxWxH



3 Tape Dispenser



2 S.S. Big Packing and Preparation Table- Two way working, with 2 Overhead Shelves

Size :- 2000 x 1400 x 1500 mm LxWxH



4 Labeller :



5 Manual Continuous Band Sealer Machine



6 S.S. Instrument Trolley with 1 Bottom Shelf and 3 sided Railing
Size:- 850 x 600 x 900 mm LxWxH



7 S.S. ETO Plastic Packing Rolls Heavy duty Holder Table
Size:- 1000 x 650 x 900 + 450 mm LxWxH



8 S.S. Wrapping Paper Dispensing Trolley

Size:- 1350 x 600 x 845 mm LxWxH

9 S.S. Working Table

Size :- 1500 x 600 x 900 mm LxWxH



10 S.S. Free standing Open Storage Rack with 5 Plain Shelves

Size :- 1200 x 450 x 1800 mm LxWxH



09898622773



www.sunmedicares.com



contact@sunmedicares.com



B - 5 "Kailas", B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.





11 S.S. Wire Mesh Instrument Sterilization Trays - Dye punched

A.) Medium Size :- 340 x 250 x 70 mm LxWxH

B). Big Size :- 450 x 250 x 70 mm LxWxH

12 S.S. Storage Cupboard (Floor mounted) with inner Storage Shelves and with Lockable Transparent Doors

Size :- 900 x 440 x 1800 mm LxWxH





FOR LINEN INSPECTION AREA



1 S.S. Big Linen Inspection Table
Size :- 1800 x 1200 x 900 mm LxWxH

2 S.S. Free standing Open Storage Rack with 5 Plain Shelves
Size :- 1200 x 450 x 1800 mm LxWxH





FOR STERILE STORAGE AREA



- 1 S.S. Wire Mesh Open Storage Rack with 5 Heavy duty Wire Mesh Shelves**
Size :- 1200 x 450 x 1800 mm LxWxH



- 2 S.S. Free standing Basket Rack for holding 16 Baskets**
Size :- 1230 x 450 x 1800 mm LxWxH



- 3 S.S. SPRI Modular Stackable Sterilizing Wire rod Autoclave Basket**
A. Size:- 585 x 395 x 195 mm LxWxH
B. Size:- 585 x 395 x 95 mm LxWxH



4 S.S. Basket Trolley

Size:- 750 x 500 x 150 mm LxWxH



5 S.S. Sterile Issue Pass Box - Double Door with Manual Mechanical Interlocking and with Transparent Vision Panel

External Size :- 600 x 650 x 600 mm LxWxH



6 S.S. Sterile Goods Issue Work Table

Size :- 1200 x 600 x 850 mm LxWxH



7 S.S. Sterile Air Tight Closed Transportation Distribution Trolley

Size :- 950 x 750 x 1150 mm LxWxH

09898622773

www.sunmedicares.com

contact@sunmedicares.com

B - 5 "Kailas", B / h : A G Teachers ,
Commerce Cross Road , Iswarbhuwan,
Navrangpura ; 380009 , Gujarat.



OTHER / MISCELLANEOUS FURNITURE ITEMS

1 M.S. Staff Change Lockers for Male and Female



2 Hanging Cupboard for Hospital Cloths:





3 S.S. Round Revolving Stool with adjustable Height arrangement (without wheels & Back Rest Support)



4 S.S. Shoe Rack



5 Linen Distribution & storage Trolley



6 SS Laundry Trolley



**7 SS Dirty Utility SLUICE
HOPPER Sink.**



**8 7.SS Bed PAN Racks
(Dirty Utility)**



9 SS Sink for CSSD / Wash Area Sink (Open Sink)



10 SS Sink for CSSD / Wash Area Sink (close sink)



11 Scrub Station for OT Area